



Aeon7200


Anesthesia Machine

CE₀₁₂₃

AEOMED
Reliable Quality Thoughtful Service

- CE certified meets international safety standards and requirements
- 7" LCD screen with key operation is simple and quick, easy to use
- The built-in patented sensor with maintenance-free prevents moisture condensation, preventing damage and ensuring accurate monitoring
- The design of the rising bellows is easy to observe its running state, simple structure, safe and stable
- The divided respiratory system design makes daily maintenance more convenient
- Alloy APL valve with fast-release function reduces resuscitation time of postoperative, to avoid repeated operations
- Quick switch between manual/mechanical mode with one key, quick response, ensure patient safety

Technical Specifications

Physical & Environmental Specifications		Ventilator Specifications	
• Dimensions (H×L×W):	1386×839×628mm	• Patient type:	Adult, Pediatric
• Weight:	85kg (Basic unit)	• Ventilation modes:	IPPV, Manual
• Storage temperature:	-20~55°C	• Setting:	Touch key
• Working temperature:	10~40°C	• TV:	0~1600ml
• Storage humidity:	≤93%	• Frequency:	4~100bpm
Electrical Specifications		• I:E:	4:1~1:8
• Mains:	AC, 100~240V, 50/60Hz	• PEEP:	0~30cmH ₂ O
• Battery:	DC, 24V, ≥120min	• Inspiration pause:	OFF, 5%~60%Ti
Pneumatic Specifications		Monitoring & Alarm	
• Gas supply:	280 ~ 600kPa (40.5 ~ 87.1psi)	• Pressure:	Ppeak, Pplat, Pmean, PEEP
• Flowmeter:	Double tube in series	• Volume:	MV, Vt
• O ₂ :	0~1L/min; 1~10L/min	• Freq.	
• N ₂ O:	0~1L/min; 1~10L/min	• Gas monitoring:	FiO ₂
• Gas system:	Hypoxia protection against apnea	• Graph display:	P-t, F-t waveforms
• O ₂ flush:	25~75L/min	• Alarm:	High/Low Airway Pressure
• Breathing circuit:	Closed, Semi-closed, Semi-open		High/Low Minute Volume
• Driven mode:	Pneumatically driven and electronically controlled		High/Low FiO ₂
			Continuous high pressure
			Power failure
			O ₂ Supply failure
			Apnea
		• Alarm silence:	≤120 seconds
		Vaporizer	
		• Supports two vaporizers	
		• Optional standby vaporizer parking holder	
		• Installation mode: Selectatec® with interlock	
		• Agent: Sev. Enf. Iso. Hal.	
		• Filling type: Pour-Fill, Key-Fill, Quik-Fil®	

Remark: Above configurations include standard and option. Please check price with your Aeonmed sales representative.

AEOMED
Reliable Quality Thoughtful Service

HQ: Building 9, No.26 Outer Ring West Road, Fengtai District, Beijing 100070, China
Science Park: No. 10, Chaobai Street, Yanjiao Development Zone, Sanhe City, Hebei Province
065201, China
TEL: +86-10-58411198
Http: //www.aeonmed.com/en
FAX: +86-10-83681616-8130
E-mail: marketing@aeonmed.com

The information contained in this leaflet is correct at the date of publication. The policy of Beijing Aeonmed Co., Ltd. is one of the continued improvement to its products. Because of this policy, Aeonmed reserves the right to make any change, which may affect the information in this leaflet without giving prior notices.
AAM7200-1907